

ADDI+IVE IMPLANTS

DELIVERING UNMATCHED OPTIONS IN CERVICAL SPACER DESIGN


SureMAX-SA™ STANDALONE CERVICAL SPACER

Unmatched Options of Sizes, Lordotic Angles, Anchors and Easy to Use Instruments

7°, 10° and 14° Lordosis


Footprints of 12 x 14 mm, 14 x 16 mm, 15 x 18 mm, 15 x 20 mm

Heights of 5 mm to 12 mm



Patent Pending

SureMAX-SA™ INSTRUMENTATION




SureMAX-SA Instrument Sterilization Tray



Awl Guide Inserter will lock on to the spacer as well as provide the proper trajectory for the Straight Awl and screw placement.



Angled Shaft Instrumentation can prepare screw holes as well as drive the screw into the vertebral body when direct access is limited.



Low Profile Inserter may also be used to implant the spacer when visualization or direct access is limited.

ADDI+IVE IMPLANTS

3101 E SHEA BLVD, SUITE 122 602-795-8850 phone
PHOENIX, AZ 85028 602-595-7896 fax

Sales@AdditiveImplants.com

 ADDITIVEIMPLANTS.COM

 [ADDITIVEIMPLANTS](https://www.linkedin.com/company/additiveimplants)

