

ADDI+IVE IMPLANTS



SureMAX[®]-SA CERVICAL STANDALONE SYSTEM

SPECIFICATIONS

SureMAX[®]-SA

CERVICAL STANDALONE SYSTEM

The SureMAX[®]-SA Cervical Standalone System offers:

- + Optimized coarse open-pore structure.
- + Convex surfaces to match spinal endplates.
- + Large lateral windows to radiologically visualize bone healing.
- + Lordotic 7°, Mid-Lordotic 10°, and Hyperlordotic 14° options.
- + Four footprints: 12 x 14 mm, 14 x 16 mm, 15 x 18 mm, 15 x 20 mm.
- + Multiple fixation options including variable, fixed and helical screws.
- + Narrow low profile inserter to enhance visibility during MIS procedures.
- + Multiple straight and angled instruments for easier screw placement.

These are essential features for the next generation Cervical Interbody Spacer.
They are combined features available in only one system:

SUREMAX[®]-SA CERVICAL STANDALONE SYSTEM



SUREMAX[®]-SA SYSTEM OVERVIEW



Additive Implant's SureMAX[®]-SA Cervical Standalone System offers a wide range of implant sizes to accommodate patient anatomy as well as a comprehensive set of instrument options to support surgical technique preferences.

Spacer

- + Material: Titanium alloy (Ti6Al4V) per ASTM F3001
- + Heights from 5 mm to 12 mm (1 mm increments)
- + 4 footprints
 - Small 12 mm x 14 mm
 - Medium 14 mm x 16 mm
 - Large 15 mm x 18 mm
 - Extra Large 15 mm x 20 mm
- + 3 lordotic angles – 7, 10 and 14 deg* (*14 deg not available in 5 mm and 10 mm heights)

Screw Options

- + Material: Titanium alloy (Ti6Al4V) per ASTM F136
- + Variable and Fixed Angle Screws – 3.5 mm diameter in lengths of 14 mm, 16 mm, 18 mm, and 20 mm (T8 Hexalobe Drive Socket)
- + Variable and Fixed Angle Rescue Screws – 3.7 mm diameter in lengths of 14 mm, 16 mm and 18 mm. (T8 Hexalobe drive socket)
- + Helical Screws – 3.5 mm diameter in lengths of 14 mm, 16 mm, 18 mm and 20 mm (Left-hand Thread Drive Socket)



Screw Sizing

Spacer Footprint	15 x 20	15 x 20	15 x 18	14 x 16	12 x 14
Screw Length	20	18	18	16	14

SUREMAX®-SA SYSTEM CONTENTS

SureMAX®-SA Spacers

Depth x Width (mm)	Height (mm)	Part Numbers		
		Lordotic 7°	Mid-Lordotic 10°	Hyperlordotic 14°
Small 12 x 14	5	1022-07121405		
	6	1022-07121406	1022-10121406	1022-14121406
	7	1022-07121407	1022-10121407	1022-14121407
	8	1022-07121408	1022-10121408	1022-14121408
	9	1022-07121409	1022-10121409	1022-14121409
	10	1022-07121410	1022-10121410	1022-14121410
	11	1022-07121411	1022-10121411	1022-14121411
	12	1022-07121412	1022-10121412	1022-14121412
Medium 14 x 16	5	1022-07141605		
	6	1022-07141606	1022-10141606	1022-14141606
	7	1022-07141607	1022-10141607	1022-14141607
	8	1022-07141608	1022-10141608	1022-14141608
	9	1022-07141609	1022-10141609	1022-14141609
	10	1022-07141610	1022-10141610	1022-14141610
	11	1022-07141611	1022-10141611	1022-14141611
	12	1022-07141612	1022-10141612	1022-14141612
Large 15 x 18	5	1022-07151805		
	6	1022-07151806	1022-10151806	1022-14151806
	7	1022-07151807	1022-10151807	1022-14151807
	8	1022-07151808	1022-10151808	1022-14151808
	9	1022-07151809	1022-10151809	1022-14151809
	10	1022-07151810	1022-10151810	1022-14151810
	11	1022-07151811	1022-10151811	1022-14151811
	12	1022-07151812	1022-10151812	1022-14151812
Extra Large 15 x 20	5	1022-07152005		
	6	1022-07152006	1022-10152006	1022-14152006
	7	1022-07152007	1022-10152007	1022-14152007
	8	1022-07152008	1022-10152008	1022-14152008
	9	1022-07152009	1022-10152009	1022-14152009
	10	1022-07152010	1022-10152010	1022-14152010
	11	1022-07152011	1022-10152011	1022-14152011
	12	1022-07152012	1022-10152012	1022-14152012



SureMAX®-SA Screws

Length (mm)	Part Numbers				
	Variable Angle Bone Screw	Fixed Angle Bone Screw	Variable Angle Rescue Screw	Fixed Angle Rescue Screw	Helical Screw
14	1322-VS14	1322-FS14	1322-RVS14	1322-RFS14	1322-HB14
16	1322-VS16	1322-FS16	1322-RVS16	1322-RFS16	1322-HB16
18	1322-VS18	1322-FS18	1322-RVS18	1322-RFS18	1322-HB18
20	1322-VS20	1322-FS20			1322-HB20

SOLUTIONS BY THE PEOPLE OF ADDITIVE IMPLANTS

We are devoted to helping your patients obtain excellent surgical outcomes. We are dedicated to supporting you with novel tools, instruments and implants. We are driven by the opportunity to share our developments with the world. We are committed partners who will do everything in our power to assist you in your quest to provide the absolute best in spinal care. You can count on us to always act as ethical partners with integrity who are worthy of your trust.

ADDI+IVE IMPLANTS

3101 E SHEA BLVD, SUITE 122
PHOENIX, AZ 85028

602-795-8850 phone
602-595-7896 fax

Sales@AdditiveImplants.com

 ADDITIVEIMPLANTS.COM

 [ADDITIVEIMPLANTS](https://www.linkedin.com/company/additiveimplants)

This documentation is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part. Please refer to the package inserts for important product information, including, but not limited to, contraindications, warnings, precautions, and adverse effects.